

PRODUCT INFORMATION

PARATEX BG-355

PARATEX™ BG-355 is a high quality, high solids water based resin that is designed for easy application, quick penetration, and fast dry time. It is used as a backcoating, anticurl, hand builder, fabric stabilizer, and laminating adhesive. It dries to a firm clear thermoplastic film. It is formaldehyde free, safe, and easy to use. PARATEX™ BG-355 resists yellowing at temperatures up to 300°F.

PHYSICAL PROPERTIES

Appearance:	White semi-viscous liquid
Solids:	54 – 56%
Water Solubility:	Completely dispersible
Specific Gravity:	1.11
pH:	4.0 – 6.0
Ionic Characteristics:	Nonionic

DIRECTIONS: PARATEX™ BG-355 is used in concentrations from 2% to 50% depending on the degree of stiffness required and the desired use. When padding onto cloth as an anticurl, hand builder, or stabilizer start with a 10% solution and adjust for desired hand. Agitation may be required in the feed tank.

If kiss rolling onto the back of fabric, start with a 50% solution. If an extremely stiff backcoating is needed, use at full strength.

PARATEX™ BG-355 acts as an excellent adhesive to bond two fabrics together. Spray a light film on one of the substrates and dry thoroughly through a frame.

Clean up with water. If dried, use Rez Away™ to remove residue.

Paratex™ and Rez Away™ are trademarks of Paradigm Labs, Inc.

The data and suggested formulations in this bulletin are based on information believed to be reliable and are offered solely for evaluation, investigation and verification of the numerous factors affecting results. Paradigm Lab's products are sold with the understanding the purchasers will make their own tests to determine the suitability of these products for the particular use. We assume no liability or responsibility for any damage to person or property resulting from or incident to the use of our products.

Statements concerning the use of Paradigm Lab's products are not be construed as recommending the infringement of any patent, and no liability for infringement arising out of any such use is assumed.



PARADIGM LABS, Inc.

P.O. BOX 138, PINE GROVE, PA 17963-0138
(570) 345-2600 – FAX (570) 345-2800